Equipment	Photo	Description
Carbolite furnace		Calcined sample
Oven Memmert		Dried sample
Four decimal place weighing scales		Sample balance
Water distillation machine		Distilled water
Raw fiber extractor		Determination of raw fiber contents
Fume hood		Limit exposure to hazardous or toxic fumes, vapors or dusts
Moisture Analyzer		Analytical Moisture Analysis

Centrifuge		Separate components of a fluid complex mixture
Ultrasonic Bath		Clean many different types of surface
Calorimeter	A COLOR	Measuring the heat of chemical reactions or physical changes as well as heat capacity
pH Meters		Measuring pH
Vortex mixer VX- 100		Mix small vials of liquid
Magnetic stirrer		Stirring liquid with heating
Desiccator		Remove moisture in sample