








Equipment	Photo	Description
Oven Yamato		Dried sample
Magnetic stirrer		Stirring liquid with heating.
Euromex microscope		Observe objects
Masticator		Grinding samples
Refrigerator		Cooling storage samples at 0-10oC
Memmert oven		Keep sample at warm temperature
Two decimal place weighing scales		Sample balance

Four decimal place
weighing scales



Sample balance

Biological Safety
Cabinets Telstar



Cabinets have been especially
designed for the handling of
microbiological agents

Magnetic stirrer IKA



Mix small vials of liquid

Colony counter



Counting of bacterial & mould
colonies

Ultrasonic Bath



Clean many different types of
surface

Vortex mixer



Mix small vials of liquid

Autoclave MLS-3020



Sterile laboratory instruments and microbiological culture medium

Memmert water bath



Incubate samples in water at a constant temperature

Magnetic stirrer



Stirring liquid with heating.

pH meter



Measures the hydrogen-ion activity in water-based solutions

Infrared oven



Heating sample
